

Wee Dry[™]
LA1115



Introduction

Wee dry™ is a microprocessor-controlled, fast heating micro-tube block. It is suitable for a range of uses in life science, molecular biology, environmental and industrial laboratories and a variety of applications including incubation, enzyme reactions, and immuno-assays.

Features:

- ◆ Compact Size
- ◆ Light Weight (**800 gm**)
- ◆ Easy to use
- ◆ Convenient to carry

Technical Specification

Model	Wee Dry™
Capacity	15 x 2ml tubes
Temperature control range	RT.+5°C-100°C
Timer range	1-999min or 1-999sec
Temp. control accuracy	±0.5°C
Display accuracy	±0.1°C
Heating time (25°C to Max.)	≤15min
Temp. control accuracy	±0.3°C
Power supply	DC 24V
Dimensions (W x D x H)mm	152 × 120 × 112 mm
Weight	0.8kg



HiMedia Laboratories Pvt. Ltd.
www.himedialabs.com

- CORPORATE OFFICE -

A-516, Swastik Disha Business Park, Via Vadhani Indl Est, LBS Marg, Mumbai - 400 086, India.
Tel : 00-91-22-6147 1919 / 2500 3747 | Fax : 00-91-22-6147 1920 / 2500 5764 | Email : info@himedialabs.com