

NEW IMMUVIEW® READER



PUT IT
TO THE
TEST



HOW TO OBTAIN FAST, OBJECTIVE TEST RESULTS ...

NEW IMMUVIEW® READER – PUT IT TO THE TEST!

The new ImmuView® Reader is a portable, highly flexible, rapid and easy to use test instrument.

Not only is it straightforward to use, but it quickly transfers your results through LIS/LIMs systems.

*Build-in LIS/LIMS integration for
a smooth transfer of test results*

Simple to use interface



Objective interpretation of results **eliminates operator subjectivity**

ImmuView® Reader
and Strip Carrier

Please contact cs@ewcdiagnostics.com or call +31(0)615606050

www.ewcdiagnostics.com



WIFI/ETHERNET connection

The reader has built-in WIFI and Ethernet connectivity for extra networking and communication options.



Easy LIS/LIMS integration

The new reader has built-in LIS/LIMS integration. This way you can easily transfer the test result within your systems for the benefits of the patient.



Superior embedded software

The embedded software provides additional tools for further in-depth assay analysis and detailed test type configuration capabilities.



999

Store up to 999 tests

The reader can store up to 999 tests. You can always retrieve previous results. Before the ImmuView® Reader reach the storage limit it will let you know.

1-2-3

Easy interface

The reader has been designed with the user in mind. Just press test and you are ready to run the test.

LCD

3.4" LCD Screen

Large screen on a small reader? No problem. Not only is the LCD screen of superior quality, but you can use it wearing gloves.



Scanner and Printer

The reader can easily be configured with a scanner and printer.



Reader and strip carrier developed for SSI Diagnostica, by Axxin, a proven expert of POC IVD instrument development

Insert the ImmuView® strip in the strip carrier, and then in the reader. Press "Test" and the results are on the way. Easy.

DIAGNOSE WITH CONFIDENCE

Now you can get rapid results with minimum effort. The new ImmuView® Reader provides quick and automatic analysis of the ImmuView® *S. pneumoniae* & *Legionella* Urinary Antigen Test.


- Two rapid results in one test
- Broad *Legionella* detection
- Easy LIS integration


IMMUVIEW® DETECT MORE THAN SEROGROUP 1 OF *LEGIONELLA*

ImmuView® *Legionella* serogroups - control antigen¹

<i>Positive</i>	1 (all monoclonal subtypes)
	2, 3, 4, 5, 6, 7, 11, 12, 13, 14, 15
<i>Negative</i>	8, 9, 10

¹ Lindsay Et al: Evaluation of the rapid lateral flow ImmuView® urinary antigen test for the detection of *L. pneumophila* and *Streptococcus pneumoniae*; poster presented at ESGLI 2017





SSI Diagnostica A/S

Herredsvejen 2
3400 Hillerød
Denmark
Tel.: +45 4829 9100
www.Immuview.com

